

---

## T2D - The Gynocular Triage to Diagnose US FDA approved cloud based clinical record

The Gynocular Triage to Diagnose (T2D) Cloud based Software  
A Clinical Record Software to use together with The Gynocular™ for a structured colposcopy examination with integrated image documentation.

The T2D is a secure cloud based electronic Clinical Record that is US FDA approved and can be used wherever you are. On a smartphone, on a LapTop or stationary DeskTop at your Clinic. In the Triage to Diagnose you can fill in new patients, schedule your patients, search among the patients already archived in your T2D Software.

The Clinical record can store information within the Colposcopy management such as, patient data, demographics, history, colposcopy examination results with integrated images of examination findings, as well as all the information about the treatment and follow up. You can also send referrals from the T2D to colleagues within your organisation or answer the referrals that have been sent to you etc., all on the run if necessary.

### Benefits with The Gynocular™ Triage to Diagnose (T2D)

- Ability to capture images directly into the T2D Clinical records Software
- Can be used as a standalone system or it can be integrated into any existing EMR System
- Telemedicine capabilities great for teaching and consulting
- Save patient data safely and securely using the Cloud based software
- Access the software using any smart phone, tablet or PC
- Colposcopy management: patient demographics, history, colposcopy examination results, treatment and referrals
- No pictures are saved on the smartphone but immediately transported to cloud for safe data management
- Software available as cloud solution or one-time download instalment to your computer



---

**gynius**

www.gynius.se  
info@gynocular.com