

			
1	Product Model	MC2001 Plus	MC2002 Plus
2	Product Code	2011000	2012000
3	Monitor Type	Two way rotatable	Two way rotatable
4	Monitor Size/Ratio	3"/4:3	3"/4:3
5	View Angle	85 °	85 °
6	Monitor Rotation	Two way rotation Left-right ≥300 °;Front-back ≥120 °	Two way rotation: Left-right ≥300 °;Front-back ≥120 °
7	Camera Pixels	450000P CMOS	450000P CMOS
8	Frame Rates	30 fps	30 fps
9	Field of View(FOV)	65±10 °	65±10 °
10	Depth of Field(DOF)	10-80mm	10-80mm
11	Internal Power/Battery	18500/3.7Vdc/2040mAh Interchangeable /Rechargeable	185000/3.7Vdc/2040mAh Interchangeable /Rechargeable
12	Fully charged battery working time	≥200min	≥200min
13	Water Resistance	IPX7 for VL handle;IPX4 for the Monitor	IPX7 for VL handle, IPX4 for the Monitor
14	Internal memory	16G	16G
15	Operating temperature	0°C-50°C	0°C-50°C
16	Operating humidity	30%-80%	30%-80%
17	Storage temperature	-20°C-50°C	-20°C-50°C
18	Storage humidity	93%	93%
19	Emergency power back-up	No	No
20	Data Retrieval	Water resistance USB Port Magnetic connection	Water resistance USB Port Magnetic connection
21	Application	For pediatric patients	For adult patients
22	Specialty	Patented blade design	Patented blade design
23	Type of blades	Slimmer Disposable: Mac-00/Mac-0/Mac-1	Slimmer Disposable: Mac-2/Mac-3/Mac-3D/Mac-4
24	Product Color	Silver	Silver
25	Product Origin	Made in China	Made in China

(V1.0-20200514)