

<b>PVL-I Plus</b>			
1	<b>Product Model</b>	<b>PVL-I-01 Plus</b>	<b>PVL-I-02 Plus</b>
2	<b>Product Code</b>	2111000	2112000
3	<b>Monitor Type</b>	Two way rotatable/Completely folded	Two way rotatable/Completely folded
4	<b>Monitor Size/Ratio</b>	3"/4:3	3"/4:3
5	<b>View Angle</b>	85 °	85 °
6	<b>Monitor Rotation</b>	Two way rotation: Left-right $\geq 300^\circ$ ; Front-back $\geq 180^\circ$	Two way rotation: Left-right $\geq 300^\circ$ ; Front-back $\geq 180^\circ$
7	<b>Camera Pixels</b>	960 X 720P CMOS	960 X 720P CMOS
8	<b>Frame Rates</b>	30 fps	30 fps
9	<b>Field of View(FOV)</b>	65±10 °	65±10 °
10	<b>Depth of Field(DOF)</b>	10-80mm	10-80mm
11	<b>Internal Power/Battery</b>	18500/3.7Vdc/2040mAh/ Interchangeable /Rechargeable	18500/3.7Vdc/2040mAh/ Interchangeable /Rechargeable
12	<b>Fully charged battery working time</b>	$\geq 200$ min	$\geq 200$ min
13	<b>Water Resistance</b>	IPX7 for VL handle, IPX4 for the Monitor	IPX7 for VL handle IPX4 for the Monitor
14	<b>Internal memory</b>	16G	16G
15	<b>Operating temperature</b>	0°C-50°C	0°C-50°C
16	<b>Operating humidity</b>	30%-80%	30%-80%
17	<b>Storage temperature</b>	-20°C-50°C	-20°C-50°C
18	<b>Storage humidity</b>	93%	93%
19	<b>Emergency power back-up</b>	No	No
20	<b>Data Retrieval</b>	Water resistance USB Port Magnetic connection	Water resistance USB Port Magnetic connection
21	<b>Application</b>	For pediatric patients	For adult patients
22	<b>Specialty</b>	Pocket Portable/patented blade	Pocket Portable/patented blade
23	<b>Type of blades</b>	Slimmer Disposable: Mac-00/Mac-0	Slimmer Disposable; Mac-1/Mac-2/Mac-3 /Mac-4
24	<b>Product Color</b>	Black	Black
25	<b>Product Origin</b>	<b>Made in the Netherlands</b>	Made in the Netherlands

(V1.0-20200514)